

Video Laryngoscope UED-C series



UED-C



- Infection control
- Cost control

EUMed

EUMed Sp. z o. o.
ul. Chłopickiego 50, 04-275 Warszawa
tel.: 022 879 69 20, fax: 022 879 69 28,
e-mail: info@eumed.pl www.eumed.pl

ZHEJIANG UE MEDICAL CORP.


✉ EXPORT@UEWORLD.COM

☎ +86 571 28057363

🌐 WWW.UEWORLD.COM/EN.PHP

📍 8 YOUYI RD, BAITA ECONOMIC
DEVELOPMENT ZONE, XIANJU,
ZHEJIANG 317300, CHINA

Signature Blades Design according to airway anatomy of different ages



UED-M3-S adds on



Touch screen



Playback function



Wi-Fi ready



Auto shutdown

Compatible with UE M8 8" monitor For various application scenarios



SPECIFICATIONS	
Camera resolution	7.87lp/mm
Angle of View	60°
Illuminance	≥150lux
Tilt and Rotation	110° tilt, 270° rotate
Memory	32G



Anti-fog tech
No need to preheat

